



NOMINAL DIMENSIONS IN MM (IN.)					
SIZE	IMP. SIZE	ØA	B	ØC	ULTRAPHARMA P/N
DN12 (050)	½"	34,0 (1,339)	16,0 (0,630)	18,0 (0,709)	UP-RD-MPE-DN12
DN18 (075)	¾"	46,5 (1,831)	22,5 (0,886)	25,0 (0,984)	UP-RD-MPE-DN18
DN25 (100)	1"	62,0 (2,441)	28,0 (1,102)	34,0 (1,339)	UP-RD-MPE-DN25
DN38 (150)	1½"	75,0 (2,953)	41,0 (1,614)	48,0 (1,890)	UP-RD-MPE-DN38
DN51 (200)	2"	97,0 (3,819)	51,0 (2,008)	64,0 (2,520)	UP-RD-MPE-DN51





SECTION VIEW


DN51 REPLACEMENT DIAPHRAGM
MATERIAL ELASTOMER: EPDM
MATERIAL METAL STEM: 316SS

Complies with CFR 177.2600 and Conforms to USP Class VI Biological Reactivity Tests Section #87 and #88	 Ultrapharma _{nl}	CONFIGURATION	DESCRIPTION:				
		EPDM	DIAPHRAGM DN51 (200) ASSEMBLY				
		The information contained in this drawing is the sole property of Ultrapharma BV. Any reproduction in part or as a whole without written permission of Ultrapharma BV is prohibited.		3RD ANGLE PROJECTION 	UNITS: mm	SCALE: 1:1	SHEET SIZE A3
info@ultrapharma.com www.ultrapharma.com		De Drogmakerij 12 1851LX Heiloo		The Netherlands Tel.: +31 856 200 900			

NOMINAL DIMENSIONS IN MM (IN.)					
SIZE	IMP. SIZE	ØA	B	ØC	ULTRAPHARMA P/N
DN12 (050)	½"	34,0 (1,339)	16,0 (0,630)	18,0 (0,709)	UP-RD-RXPXC-DN12
DN18 (075)	¾"	46,5 (1,831)	22,5 (0,886)	25,0 (0,984)	UP-RD-RXPXC-DN18
DN25 (100)	1"	62,0 (2,441)	28,0 (1,102)	34,0 (1,339)	UP-RD-RXPXC-DN25
DN38 (150)	1½"	75,0 (2,953)	41,0 (1,614)	48,0 (1,890)	UP-RD-RXPXC-DN38
DN51 (200)	2"	97,0 (3,819)	51,0 (2,008)	64,0 (2,520)	UP-RD-RXPXC-DN51



SECTION VIEW



DN51 REPLACEMENT DIAPHRAGM
MATERIAL ELASTOMER: PLAT. SILICONE TRANSLUCENT
MATERIAL METAL STEM: 316SS

Complies with CFR 177.2600
and Conforms to
USP Class VI Biological
Reactivity
Tests Section #87 and #88


UP Ultrapharma

info@ultrapharma.com De Droogmakerij 12 The Netherlands
www.ultrapharma.com 1851LX Heiloo Tel.: +31 856 200 900

CONFIGURATION
PLAT. SILICONE TRANSLUCENT

DESCRIPTION:
DIAPHRAGM DN51 (200) ASSEMBLY

The information contained in this drawing is the sole
property of Ultrapharma BV. Any reproduction in
part or as a whole without written permission of
Ultrapharma BV is prohibited.





UNITS:
mm



SCALE:
1:1

SHEET SIZE
A3

NOMINAL DIMENSIONS IN MM (IN.)					
SIZE	IMP. SIZE	ØA	B	ØC	ULTRAPHARMA P/N
DN12 (050)	½"	34,0 (1,339)	16,0 (0,630)	18,0 (0,709)	UP-RD-RXPXW-DN12
DN18 (075)	¾"	46,5 (1,831)	22,5 (0,886)	25,0 (0,984)	UP-RD-RXPXW-DN18
DN25 (100)	1"	62,0 (2,441)	28,0 (1,102)	34,0 (1,339)	UP-RD-RXPXW-DN25
DN38 (150)	1½"	75,0 (2,953)	41,0 (1,614)	48,0 (1,890)	UP-RD-RXPXW-DN38
DN51 (200)	2"	97,0 (3,819)	51,0 (2,008)	64,0 (2,520)	UP-RD-RXPXW-DN51



SECTION VIEW



DN51 REPLACEMENT DIAPHRAGM
MATERIAL ELASTOMER: PLAT. SILICONE WHITE
MATERIAL METAL STEM: 316SS

Complies with CFR 177.2600
 and Conforms to
 USP Class VI Biological
 Reactivity
 Tests Section #87 and #88


UP Ultrapharma

info@ultrapharma.com De Droogmakerij 12 The Netherlands
 www.ultrapharma.com 1851LX Heiloo Tel.: +31 856 200 900

CONFIGURATION
 PLAT. SILICONE WHITE

DESCRIPTION:
 DIAPHRAGM DN51 (200) ASSEMBLY

The information contained in this drawing is the sole
 property of Ultrapharma BV. Any reproduction in
 part or as a whole without written permission of
 Ultrapharma BV is prohibited.



UNITS:
 mm

SCALE:
 1:1

SHEET SIZE
 A3